

# **Absorbance Theory**

## **- Review -**

*How to measure absorbance using a microplate reader  
and ELISA Basics*

### Goals

---

- Understand absorbance numbers

1.572

2.853

0.143

0.042

- Concepts

- Absolute vs. Relative

- Quantify vs. Qualify

- Understand why wavelength selection is important

- How BioTek verifies the reader is working correctly

- ELISA Basics

### First two concepts

---

- Absolute vs. Relative
  - Absolute = Absorbance
  - Relative = Fluorescence, Luminescence
  
- Quantify vs. Qualify
  - Quantify = Amount
  - Qualify = Good or Bad (or in between)


### What is Optical Absorption?

---

- Optical absorption is a process where **light energy** is **transformed** by a molecule into another form of energy, usually **heat**.
- To understand Absorbance, it's helpful to first understand transmittance.

## What is transmittance?

---





Path length

$$T = \frac{I}{I_0}$$

It's a Ratio


You need to measure **both**  $I$  and  $I_0$  to calculate the %T

## What is Absorbance (AU)?



**Absorbance units (AU) are also called optical density (OD)**

% Transmittance	Absorbance
100	0.000
90	0.046
80	0.097
70	0.155
60	0.222
50	0.301
40	0.398
30	0.523
20	0.699
10	1.000
5	1.301
2	1.699
1	2.000
0.5	2.301
0.1	3.000
0.01	4.000
0.001	5.000




## Light as a wave

# Is color important?


- **Wavelength**

- **Amplitude**






### Electromagnetic Radiation




### Electromagnetic Radiation



Color	Wavelength (nm) Approximate
Ultraviolet	001-380
Violet	380-430
Blue	430-500
Cyan	500-520
Green	520-565
Yellow	565-590
Orange	590-625
Red	625-740
Infrared	750-1,000,000

### An object's color

- When white light falls on an object, the object **absorbs** certain waves and reflects others, this determines the object's color.



### Importance of color in the Assay

---

•What color is the liquid in the microplate?

•What wavelength are you using to read the sample?

•Why?

Color	Wavelength (nm) Approximate
Ultraviolet	001-380
Violet	380-430
Blue	430-500
Cyan	500-520
Green	520-565
Yellow	565-590
Orange	590-625
Red	625-740
Infrared	750-1,000,000

### Components of the absorbance instrument

---

- |                                |                   |
|--------------------------------|-------------------|
| 1. Light source                | Tungsten Halogen  |
| 2. Wavelength selection device | Band Pass Filters |
| 3. Detector                    | Photodiode        |
| 4. Method to acquire $I_0$     | Computer          |

## Qualification

---

- Absorbance Test Plate
  - Solid state neutral density glass
- Liquid Tests
  - Color solution

Farayano



### Qualification

---

#### Absorbance Test Plate

- Alignment
- Accuracy
- Repeatability
- Turnaround
- Linearity
- Wavelength Accuracy (Not use with the ELx808)

## Qualification


---

### Liquid Tests

- Liquid testing differs from testing with the Absorbance Test Plate because liquid in the wells has a **meniscus**.
- The optics characteristics may differ in these two methods.



### Liquid Tests



**The curved meniscus acts like a lens.**

Microplate (Flat Bottom)

Fluid used is DI water from 50ul to 300ul

**To reduce this problem, we recommend shaking the plate for 5 minutes or letting it sit 20 minutes before reading the plate.**

### Liquid Tests

---

- Liquid Test 1
  - Alignment, Accuracy, Repeatability, and channel-to-channel variability
- Liquid Test 2 (Not needed if using Absorbance Test Plate)
  - Linearity, Repeatability, and Alignment
  - Uses a series of solutions of varying absorbencies
- Liquid Test 3
  - An optional test offered for sites that must have proof of linearity at wavelengths lower than those attainable with the Absorbance Test Plate
  - This test is considered optional because the reader has good “front end” linearity throughout its wavelength range

### Liquid Test 1

#### Absorbance Liquid Test 1

Mean ODs  
for five reads:

"Normal"

	1	2
A	1.988	1.023
B	2.008	1.023
C	2.006	1.031
D	2.003	1.024
E	2.010	1.034
F	2.016	1.032
G	2.016	1.029
H	2.009	1.030

"Turnaround"

	11	12
A	1.029	2.004
B	1.027	2.021
C	1.029	2.011
D	1.034	2.020
E	1.025	2.010
F	1.034	2.022
G	1.022	2.010
H	1.019	1.999

Results:

Wells for Comparison	Spec	Range		Result
		Low	High	
A1/H12	+/- 0.030	1.9581	2.0179	PASS
B1/G12	+/- 0.030	1.9779	2.0381	PASS
C1/F12	+/- 0.030	1.9759	2.0361	PASS
D1/E12	+/- 0.030	1.9730	2.0330	PASS
E1/D12	+/- 0.030	1.9799	2.0401	PASS
F1/C12	+/- 0.030	1.9858	2.0462	PASS
G1/B12	+/- 0.030	1.9858	2.0462	PASS
H1/A12	+/- 0.030	1.9789	2.0391	PASS
A2/H11	+/- 0.020	1.0028	1.0432	PASS
B2/G11	+/- 0.020	1.0028	1.0432	PASS
C2/F11	+/- 0.020	1.0107	1.0513	PASS
D2/E11	+/- 0.020	1.0038	1.0442	PASS
E2/D11	+/- 0.020	1.0137	1.0543	PASS
F2/C11	+/- 0.020	1.0117	1.0523	PASS
G2/B11	+/- 0.020	1.0087	1.0493	PASS
H2/A11	+/- 0.020	1.0097	1.0503	PASS


Accuracy specification:

± 1.0% ± 0.010 Abs from 0 to 2.000 Abs

### Liquid Test 3

Average Mean for each dye concentration ('Y'):

	X	Y
Wells A5-H6 Low-level (50%):	0.50	0.511
Wells A3-H4 Mid-level (75%):	0.75	0.704
Wells A1-H2 High-level (100%):	1.00	0.910



\* Also checks  
Standard  
Deviation

## **ELISA**

---

- **Definitions**
- **Enzyme-Linked Immunosorbent Assay**
- **ELISA Applications**
- **[Sumanas ELISA Video](#)**

## ELISA Definitions

---

### Antibody

Antibodies are a type of [protein](#). They are produced by the immune system in response to foreign substances that may be a threat to the body -- such as chemicals, virus particles, spores, or bacterial toxins. These foreign substances are called antigens.

Each type of antibody is unique and defends the body against one specific type of antigen.

### Antigen

An antigen is any substance that causes your immune system to produce [antibodies](#) against it. The antigen may be a foreign substance from the environment (such as chemicals, bacteria, viruses, or pollen) or formed within the body (such as bacterial [toxins](#) or tissue cells).

### Enzyme

Enzymes are complex proteins that cause a specific chemical change in other substances without being changed themselves. Enzymes can be found in every organ of the body. For example, they can change [starches](#), [proteins](#), and [sugars](#) into substances the body can digest. Blood clotting is another example of enzymes at work.

Enzymes exist in the mouth (saliva), stomach (gastric juice), and intestines (pancreatic juice, intestinal juice, and intestinal [mucosa](#)).

## ELISA

---

- **Enzyme-Linked Immunosorbent Assay**

- Alternate Names : EIA, Enzyme-Linked Immunoassay

- This versatile test is widely used in the clinical laboratory

- It allows your health care provider to:

- test your blood with an antigen (e.g., virus or bacteria) to see if your immune system recognizes it as something it has seen before, or

- test your blood with an antibody to see if a particular substance like a hormone (an antigen) is present in your system.

## ELISA

---

### ELISA Applications

- HIV
- Lyme disease
- Cell death detection
- Cell Proliferation
- Blood Donor Screening
- Syphilis
- Hepatitis
- Bird Flu
- H1N1
- hCG Pregnancy Test (Human Chorionic Gonadotropin)
- Anabolic Steroids
- Narcotics
- LSD
- Drug Facilitated Sexual Assault

The list goes on and on....



## **ELISA**

---

**[Sumanas ELISA Video](#)**

Farayand



## Questions

---

Questions?

Thank You

Farooq and

